

Forehead Infrared Thermometer



Model No.: FC-IR209



Three Color
Indicator Light



Forehead
Measure



One Second
Reading



Mute Mode



Fever Warning



White BackLight



35 Memory
Function



Clinically Tested



Object/Room
Measure



Infrared
Technology

Features

- Instant reading; Quick measurement, get reading in 1 second.
- Adopt the latest algorithm with optimised temperature calibration function.
- Accurate and reliable; $\pm 0.2^{\circ}\text{C}$ High accuracy.
- Smart and Intelligent; One-button easy operation.
- Tri-color indicator light on measure button, display green, orange & red with fever alarm function.
- Memory function. Memory recall of 35 reading.
- Multi function design, can measure forehead/ear, room, object temperature.
- Switching between mute and un-mute.
- Switching between Fahrenheit and Celsius.
- Automatic shutdown and power saving.
- Optional function: Bluetooth connection.

Specifications

Model No.	FC-IR209
Measurement range	Forehead: 32.0°C – 42.9°C 89.6°F – 109.2°F Object: 0°C – 100°C (32°F – 212°F)
Measuring accuracy	Forehead mode: $\pm 0.2^{\circ}\text{C}$ / $\pm 0.4^{\circ}\text{F}$ Object mode: $\pm 1.0^{\circ}\text{C}$ / 1.8°F
Display	LED backlight
Automatic shutdown	10 seconds automatic shutdown function
Memory	Memory 35 groups of measured temperature.
Indicator and Fever alarm	Normal green: $32^{\circ}\text{C} \leq T \leq 37.3^{\circ}\text{C}$ (89.6°F – 109.2°F) Low fever: $37.4^{\circ}\text{C} \leq T \leq 37.9^{\circ}\text{C}$ (99.3°F – 100.3°F) High fever: $38^{\circ}\text{C} \leq T \leq 42.5^{\circ}\text{C}$ (100.4°F – 108.6°F)
Operational conditions	Temperature: 10°C – 40°C (50°F – 104°F) Humidity: 15–95%RH, non-condensing Atmospheric pressure: 86–106 kPa
Battery	2*1.5V(AAA)
Weight	74g (without battery)
Product dimension	148.6*35.5*42mm
Packing size	46*46*168 mm
Full set weight	121g(including products/bags/instructions/battery)