



Ear Temperature measuring


Infrared thermometer FC-IR202



- Product dimension: 36*42*153.5mm
- Packing size: 46*46*168mm
- Full set weight: 115g
- Weight of thermometer: 66.8g(Without battery)/81.4g
- Quantity per carton: 100 Pieces
- N.W/Carton: 12.5kg
- G.W/Carton: 14kg
- Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch 2*Battery, 1*Color box




Infrared thermometer FC-IR1030



- Product dimension: 45*53*120mm
- Packing size: 46*55*127mm
- Full set weight: 106.6g
- Weight of thermometer: 60.8g(Without battery)/75.4g
- Quantity per carton: 100 Pieces
- N.W/Carton: 11.5kg
- G.W/Carton: 13kg
- Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch 2*Battery, 1*Color box



Blood pressure monitor FC-BP113



- Large LED display size: 120×100mm
- Accuracy: Pressure within ±3mmHg(±0.4kpa), pulse within ±5%.
- Optional: Voice broadcast; Bluetooth; LED with backlight.
- Power supply: 4×AA alkaline batteries (DC6V); USB Type C (DC5V1A)
- Product weight: 335g
- Product dimension: 112×155×79mm
- Packing size: 118×110×160mm
- Cuff size: 22-42cm

Blood pressure monitor FC-BP121



Infrared thermometer FC-IR109



- Product dimension: 35*69.4*145.5mm
- Packing size: 45*80*165mm
- Full set weight: 129.4g
- Weight of thermometer: 70g(Without battery)/85.6g
- Quantity per carton: 100 Pieces
- N.W/Carton: 14.5kg
- G.W/Carton: 16kg
- Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch 2*Battery, 1*Color box



Infrared thermometer FC-IR207



- Product dimension: 157*35*36mm
- Packing size: 46*46*168mm
- Full set weight: 136g
- Weight of thermometer: 69g(Without battery)/85g
- Quantity per carton: 100 Pieces
- N.W/Carton: 9kg
- G.W/Carton: 10.5kg
- Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch 2*Battery, 1*Color box



Blood pressure monitor FC-BP100



- Large LED display size: 120×100mm
- Accuracy: Pressure within ±3mmHg(±0.4kpa), pulse within ±5%.
- Optional: Voice broadcast; LED with backlight.
- Power supply: 4×AA alkaline batteries (DC6V); USB Type C (DC5V1A)
- Product weight: 298g
- Product dimension: 109×144×64mm
- Packing size: 110×95×160mm
- Cuff size: 22-42cm

Blood pressure monitor FC-BP130



- Large LED display size: 120×100mm
- Accuracy: Pressure within ±3mmHg(±0.4kpa), pulse within ±5%.
- Optional: Voice broadcast; Bluetooth; LED with backlight.
- Power supply: 4×AA alkaline batteries (DC6V); USB Type C (DC5V1A)
- Product weight: 248g
- Product dimension: 132×110.6×66mm
- Packing size: 120×98×159mm
- Cuff size: 22-42cm

Infrared thermometer FC-IR100



- Product dimension: 39*27.3*160mm
- Packing size: 46*46*168mm
- Full set weight: 113g
- Weight of thermometer: 66g(Without battery)/81g
- Quantity per carton: 100 Pieces
- N.W/Carton: 12.5kg
- G.W/Carton: 14kg
- Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch 2*Battery, 1*Color box




Infrared thermometer FC-IR1010



- Product dimension: 36*42*153.5mm
- Packing size: 46*46*168mm
- Full set weight: 115g
- Weight of thermometer: 66.8g(Without battery)/81.4g
- Quantity per carton: 100 Pieces
- N.W/Carton: 12.5kg
- G.W/Carton: 14kg
- Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch 2*Battery, 1*Color box




Blood pressure monitor FC-BP103




- Large LED display size: 120×100mm
- Accuracy: Pressure within ±3mmHg(±0.4kpa), pulse within ±5%.
- Optional: Voice broadcast; LED with backlight.
- Power supply: 4×AA alkaline batteries (DC6V); USB Type C (DC5V1A)
- Product weight: 298g
- Product dimension: 109×144×64mm
- Packing size: 110×95×160mm
- Cuff size: 22-42cm

Blood pressure monitor FC-BP131



- Large LED display size: 120×100mm
- Accuracy: Pressure within ±3mmHg(±0.4kpa), pulse within ±5%.
- Optional: Voice broadcast; Bluetooth; LED with backlight.
- Power supply: 4×AA alkaline batteries (DC6V); USB Type C (DC5V1A)
- Product weight: 249g
- Product dimension: 132×110.6×66mm
- Packing size: 120×98×159mm
- Cuff size: 22-42cm



- Large LED display size: 120×100mm
- Accuracy: Pressure within ±3mmHg(±0.4kpa), pulse within ±5%.
- Optional: Voice broadcast; Bluetooth; LED with backlight.
- Power supply: 4×AA alkaline batteries (DC6V); USB Type C (DC5V1A)
- Product weight: 335g
- Product dimension: 129×96×53mm
- Packing size: 118×110×160mm
- Cuff size: 22-42cm

Certificates

TGA CE0123 TUV ISO13485 RoHS

FDA CFDA BSCI FCC FSC



CFDA No.: 20202071679

FDA 510K No.: K181215



深圳市福瑞诺科技有限公司

Shenzhen Finicare Co., Ltd.



OUR VISION

Lead the industry, to be a international leader in manufacturing and services of high end quality of medical devices.

Finicare, a step towards health care



Web



Wechat

SHENZHEN FINICARE CO., LTD

sales@finicare.com

info@finicare.com

201, No. 50, the 3rd Industrial Park Houting Community, Shajing Street Bao'an District 518104 Shenzhen

https://finicare.com/

+86 755 2301 3503

+86 135 3051 9669

+86 177 2759 8669

About FINICARE

Shenzhen Finicare, a professional manufacturer of medical device with FDA/CE accreditation, headquarter in Shenzhen, China.

Finicare is a world's leading supplier and manufacturer of medical device with a R&D team of more than 14 years of experience featuring infrared technique. Finicare attaches great importance to independent research, boasting 35 patent technologies and software copyrights.

Finicare have developed multiple product lines including Infrared thermometer, Blood pressure monitor, Digital thermometer, Mesh Nebulizer and Oximeter, etc. With implementation of ISO13485 and establishment of ERP system, we have formulated a complete set of management system ranging from research, development, production to quality control and secure shipment. Our products have passed certifications of FDA, CE(TUV), CFDA, FSC, FCC, ROHS, BSCI and TGA, which is the best guarantee of quality of our product.

Currently, our products are exported to more than 80 countries and regions including America, Europe, Australia, Asia, Africa, etc. We have established exclusive agency partnerships with several countries such as Netherlands, Saudi Arabia, Belgium and developed stable markets in other places. We have achieved considerable development in terms of our own brands and OEM production. We are persistent in our pursuit of improving our current 'Standard of care'-constantly researching, testing and reinventing conventional processes in order to make the world a better place with our healthcare technology. Our company's capabilities of R&D, production, and quality assurance have been recognized by more than 650,000 customers worldwide.

Finicare will unremittingly apply the most advanced technology to R&D of medical device for meeting high demand we are persisting in and making our world-class product.

Finicare covers an area of 5,000 square meters and has more than 200 employees. Our production capability is 500,000 pieces per month. Finicare firmly uphold the idea that quality is the soul of product, which should be strictly implemented throughout the whole process ranging from implemetation of product selection standard, product manufacturing, guarantee of service detail, commitment of maximum value to customers. Finicare is renowned for its high quality and reliable product characterized by best performance.

Our mission:
We are committed in persistent pursuit of providing the best medical device and service for our customers worldwide.

Our vision:
Becomes an intermation leader in the healthcare industry by providing sophisticated and high quality medical device and service.

