



深圳市福瑞諾科技有限公司

PRODUCT CATALOG

Finicare, a step towards health care

Contents

| | | |
|-----------|---|--------------|
| 1) | Company introduction | 02 |
| 2) | Certificates | 03 |
| 3) | Digital Blood Pressure Monitor | 04-20 |
| 4) | Infrared thermometer | 21-38 |

Introduction

Shenzhen Finicare, a professional manufacturer of medical device with FDA/CE accreditation, headquarter in Shenzhen, China.

Finicare is a world' s leading supplier and manufacturer of medical device with a R&D team of more than 14 years of experience featuring infrared technique. Finicare attaches great importance to independent research, boasting 35 patent technologies and software copyrights.

Finicare have developed multiple product lines including Infrared thermometer, Blood pressure monitor, Digital thermometer , Mesh Nebulizer and Oximeter,etc. With implementation of ISO13485 and establishment of ERP system, we have formulated a complete set of management system ranging from research, development, production to quality control and secure shipment. Our products have passed certifications of FDA, CE(TUV), CFDA, FSC, FCC, ROHS, BSCI and TGA, which is the best guarantee of quality of our product.

Currently, our products are exported to more than 80 countries and regions including America, Europe, Australia, Asia, Africa,etc. We have established exclusive agency partnerships with several countries such as Netherlands, Saudi Arabia, Belgium and developed stable markets in other places. We have achieved considerable development in terms of our own brands and OEM production.

We are persistent in our pursuit of improving our current 'Standard of care' -constantly researching, testing and reinventing conventional processes in order to make the world a better place with our healthcare technology. Our company' s capabilities of R&D, production, and quality assurance have been recognized by more than 650,000 customers worldwide.

Finicare will unremittingly apply the most advanced technology to R&D of medical device for meeting high demand we are persisting in and making our world-class product.

ADVANTAGES:

Finicare covers an area of 5,000 square meters and has more than 200 employees. Our production capability is 500,000 pieces per month. Finicare firmly uphold the idea that quality is the soul of product, which should be strictly implemented throughout the whole process ranging from implemetation of product selection standard, product manufacturing, guarantee of service detail, commitment of maximum value to customers. Finicare is renowned for its high quality and reliable product characterized by best performance.

Our mission:

We are committedly in persistent pursuit of providing the best medical device and service for our customers worldwide.

Our vision:

Becomes an intarnation leader in the healthcare industry by providing sophisticated and high quality medical device and service.



Certificates



FDA 510K



CE



TGA



ISO 13485



BSCI



FCC



ROHS



CFDA

Blood Pressure Monitor



FC-BP100LED

FC-BP100LED



Easy 1-touch Operation




WHO Indicator



Automatically Power Off



Noble Black Design



Intelligent Measurement



Arrhythmia Detection



Body Movement
Detection



LED Display



2*120 Groups Memory
Function



Type-C Connection



Self-check Cuff
Positioning



Colorful Display

Model No: FC-BP100LED

Large LED display size: 120mm×100mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; LED with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 298g

Product dimension: 109×144×64mm

Packing size: 110×95×160mm

Cuff size: 22-42cm



Blood Pressure Monitor



FC-BP103

FC-BP103



Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Noble Black Design



Intelligent Measurement



Arrhythmia Detection



Body Movement
Detection



LED Display



2*120 Groups Memory
Function



Type-C Connection



Self-check Cuff
Positioning



Colorful Display

Model No: FC-BP103

Large LED display size: 120mm×100mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; LED with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 298g

Product dimension: 109×144×64mm

Packing size: 110×95×160mm


Cuff size: 22-42cm

Blood Pressure Monitor



FC-BP113

FC-BP113



Model No: FC-BP113

Large LED display size: 120mm×100mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; Bluetooth; LED with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 335g

Product dimension: 112×155×79mm

Packing size: 118×110×160mm

Cuff size: 22-42cm



Finicare, a step towards health care

Blood Pressure Monitor



FC-BP121

FC-BP121



Easy 1-touch Operation




WHO Indicator



Automatically Power Off



Noble Black Design



Intelligent Measurement



Arrhythmia Detection



Body Movement
Detection



LED Display



2*120 Groups Memory
Function



Type-C Connection



Self-check Cuff
Positioning



Colorful Display

Model No: FC-BP121

Large LED display size: 80mm×75mm

Accuracy: Pressure within $\pm 3\text{mmHg}$ ($\pm 0.4\text{kpa}$)

Pulse within $\pm 5\%$

Optional: Voice broadcast; Bluetooth; LED with backlight

Power supply: 4×AAA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 335g

Product dimension: 129×96×53mm

Packing size: 118×110×160mm

Cuff size: 22-42cm



Blood Pressure Monitor



FC-BP130

FC-BP130



Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Noble Black Design



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



LED Display



2*120 Groups Memory Function



Type-C Connection



Self-check Cuff Positioning



Colorful Display

Model No: FC-BP130

Large LED display size: 78mm×78mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; Bluetooth; LED with backlight

Power supply: 4×AAA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 247.8g

Product dimension: 132 x 110.6 x 66mm

Packing size: 120 x 98 x 159mm

Cuff size: 22-42cm

Blood Pressure Monitor



FC-BP131

FC-BP131



Easy 1-touch Operation




WHO Indicator



Automatically Power Off



Noble Black Design



Intelligent Measurement



Arrhythmia Detection



Body Movement
Detection



LED Display



2*120 Groups Memory
Function



Type-C Connection



Self-check Cuff
Positioning



Colorful Display

Model No: FC-BP131

Large LED display size: 78mm×78mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; Bluetooth; LED with backlight

Power supply: 4×AAA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 249.2g

Product dimension: 132 x 110.6 x 66mm

Packing size: 120 x 98 x 159mm

Cuff size: 22-42cm

Blood Pressure Monitor CE FDA



FC-BP120

FC-BP120




Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*99 Groups Memory
Function



Type-C Connection



Self-check Cuff
Positioning

Model No: FC-BP120

Large LCD display size: 60.8×75.3mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; LCD with backlight

Power supply: 4×AAA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 240g

Product dimension: 129×96×53mm

Packing size: 118×110×160mm

Cuff size: 22-42cm



Finicare, a step towards health care

Blood Pressure Monitor



FC-BP100

FC-BP100




Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory Function



Type-C Connection



Self-check Cuff Positioning

Model No: FC-BP100

Large LCD display size: 60.8×75.3mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; Bluetooth; LCD with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 299g

Product dimension: 109×144×64mm

Packing size: 110×95×160mm

Cuff size: 22-42cm



Blood Pressure Monitor



FC-BP101

FC-BP101



Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory
Function



Type-C Connection



Self-check Cuff
Positioning

Model No: FC-BP101

Large LCD display size: 60.8×75.3mm

Accuracy: Pressure within $\pm 3\text{mmHg}$ ($\pm 0.4\text{kpa}$)

Pulse within $\pm 5\%$

Optional: Voice broadcast; LCD with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 276g

Product dimension: 109×144×64mm

Packing size: 110×95×160mm

Cuff size: 22-42cm



Blood Pressure Monitor



FC-BP102

FC-BP102



Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory
Function



Type-C Connection



Self-check Cuff
Positioning

Model No: FC-BP102

Large LCD display size: 60.8×75.3mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; LCD with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 265g

Product dimension: 109×144×64mm

Packing size: 110×95×160mm

Cuff size: 22-42cm



Blood Pressure Monitor



FC-BP110

FC-BP110




Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory Function



Type-C Connection



Self-check Cuff Positioning

Model No: FC-BP110

Large LCD display size: 60.8×75.3mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; Bluetooth; LCD with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 312g

Product dimension: 105×150×76mm

Packing size: 110×95×160mm

Cuff size: 22-42cm



Finicare, a step towards health care

Blood Pressure Monitor



FC-BP111

FC-BP111




Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory
Function



Type-C Connection



Self-check Cuff
Positioning

Model No: FC-BP111

Large LCD display size: 60.8×75.3mm

Accuracy: Pressure within $\pm 3\text{mmHg}$ ($\pm 0.4\text{kpa}$)

Pulse within $\pm 5\%$

Optional: Voice broadcast; Bluetooth; LCD with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 365g

Product dimension: 112×155×79mm

Packing size: 118×110×160mm

Cuff size: 22-42cm

Blood Pressure Monitor

FDA



FC-BP112

FC-BP112




Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory Function



Type-C Connection



Self-check Cuff Positioning

Model No: FC-BP112

Large LCD display size: 60.8×75.3mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Voice broadcast; LCD with backlight

Power supply: 4×AA alkaline batteries (DC6V)

USB Type C (DC5V1A)

Product weight: 325g

Product dimension: 112×155×78mm

Packing size: 118×110×160mm

Cuff size: 22-42cm



Blood Pressure Monitor



FC-BP200

FC-BP200



Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory
Function



Optional: Voice Broadcast



Optional: Bluetooth

Model No: FC-BP200

Large LCD display size: 47×51mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Bluetooth; Voice Broadcast

Power supply: 2 AAA alkaline batteries

(not included), DC3V

Product weight: 95.4g

Product dimension: 80.4×65.5×24mm

Packing size: /

Cuff size: 13.5-21.5cm

Blood Pressure Monitor



FC-BP210

FC-BP210




Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory
Function



Optional: Voice Broadcast



Optional: Bluetooth

Model No: FC-BP210

Large LCD display size: 41×43mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Bluetooth; Voice Broadcast

Power supply: 2 AAA alkaline batteries

(not included), DC3V

Product weight: 95.4g

Product dimension: 80.4×65.5×24mm

Packing size: /

Cuff size: 13.5-21.5cm



Finicare, a step towards health care

Blood Pressure Monitor



FC-BP220

FC-BP220




Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory
Function



Optional: Voice Broadcast



Optional: Bluetooth

Model No: FC-BP220

Large LCD display size: 76×58mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Bluetooth; Voice Broadcast

Power supply: 2 AAA alkaline batteries

(not included), DC3V

Product weight: 95.4g

Product dimension: 96×63×20mm

Packing size: /

Cuff size: 13.5-21.5cm



Blood Pressure Monitor



FC-BP221

FC-BP221




Easy 1-touch Operation



WHO Indicator



Automatically Power Off



Intelligent Measurement



Arrhythmia Detection



Body Movement Detection



2*120 Groups Memory Function



Optional: Voice Broadcast



Optional: Bluetooth

Model No: FC-BP221

Large LCD display size: 76×58mm

Accuracy: Pressure within $\pm 3\text{mmHg}(\pm 0.4\text{kpa})$

Pulse within $\pm 5\%$

Optional: Bluetooth; Voice Broadcast

Power supply: 2 AAA alkaline batteries

(not included), DC3V

Product weight: 95.4g

Product dimension: 96×63×20mm

Packing size: /

Cuff size: 13.5-21.5cm



Infrared Thermometer CE FDA



FC-IR100

FC-IR100



Accu-sensor
Technology



Ear/Forehead
Measure



Fever Alarm



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Tri-color Backlit

Model No: FC-IR100

Product dimension: 39*27.3*160mm

Packing size: 46*46*168mm

Full set weight: 113g

Weight: 66g(Without battery)/81g

Quantity per carton: 100 Pieces

N.W/Carton: 12.5kg

G.W/Carton: 14kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch,
2*Batteries, 1*Color box



Infrared Thermometer



FC-IR209

FC-IR209



Accu-sensor
Technology



Forehead
Measure



Fever Alarm



One Second
Reading



Extra-large LED
Display



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Light Spot
Positioning



Button Aperture
Temperature
Smart Indicator

Model No: FC-IR209

Product dimension: 148.6*35.5*42mm

Packing size: 46*46*168mm

Full set weight: 121g

Weight: 74g(Without battery)/93g

Quantity per carton: 100 Pieces

N.W/Carton: 12.5kg

G.W/Carton: 14kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch,
2*Batteries, 1*Color box



Infrared Thermometer

CE FDA



FC-IR100LED

FC-IR100LED



Accu-sensor
Technology



Ear/Forehead
Measure



Fever Alarm



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Tri-color Backlit

Model No: FC-IR100LED

Product dimension: 39*27.3*160mm

Packing size: 46*46*168mm

Full set weight: 113g

Weight: 66g(Without battery)/81g

Quantity per carton: 100 Pieces

N.W/Carton: 12.5kg

G.W/Carton: 14kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch,
2*Batteries, 1*Color box



Infrared Thermometer CE FDA



FC-IR300

FC-IR300



Accu-sensor
Technology



Ear/Forehead
Measure



Baby Mode



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm

Model No: FC-IR300

Product dimension: 45*53*120mm

Packing size: 46*55*127mm

Full set weight: 106.6g

Weight of thermometer: 60.8g(Without battery)/75.4g

Quantity per carton: 100 Pieces

N.W/Carton: 11.5kg

G.W/Carton: 13kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box

Infrared Thermometer

CE FDA



FC-IR400

FC-IR400



Accu-sensor
Technology



Ear/Forehead
Measure



Fever Alarm



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Tri-color Backlit

Model No: FC-IR400

Product dimension: 43*68*161mm

Packing size: 48*68*170mm

Full set weight: 142.4g

Weight of thermometer: 85.8g(Without battery)/119.8g

Quantity per carton: 100 Pieces

N.W/Carton: 14.5kg

G.W/Carton: 16kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box



Infrared Thermometer

CE FDA



FC-IR1030

FC-IR1030



Accu-sensor
Technology



Ear/Forehead
Measure



Baby Mode



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm

Model No: FC-IR1030

Product dimension: 44*43*159mm

Packing size: 46*46*168mm

Full set weight: 116g

Weight of thermometer: 68.8g(Without battery)/83.8g

Quantity per carton: 100 Pieces

N.W/Carton: 12.8kg

G.W/Carton: 14.3kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box



Infrared Thermometer



FC-IR1030LED

FC-IR1030LED



Accu-sensor
Technology



Ear/Forehead
Measure



Baby Mode



One Second
Reading



Three Color
Indicator Light



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



LED display

Model No: FC-IR1030LED

Product dimension: 44*43*159mm

Packing size: 46*46*168mm

Full set weight: 116g

Weight of thermometer: 68.8g(Without battery)/83.8g

Quantity per carton: 100 Pieces

N.W/Carton: 12.8kg

G.W/Carton: 14.3kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box

Infrared Thermometer

CE FDA



FC-IR3010

FC-IR3010



Accu-sensor
Technology



Ear/Forehead
Measure



One Second
Reading



Fever Alarm



35 Groups Memory
Function



Object/Room
Measure



°F/°C
Switchable



Tri-color Backlit

Model No: FC-IR3010

Product dimension: 36.7*30.5*119mm

Packing size: 44*40*125mm

Full set weight: 77g

Weight of thermometer: 46.6g(Without battery)/49.8g

Quantity per carton: 100 Pieces

N.W/Carton: 9kg

G.W/Carton: 10.5kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
1*Lithium Battery, 1*Color box



Infrared Thermometer



Package includes Ear
Thermometer Probe Covers

Ear Thermometer Probe Covers Test



FC-IR109

FC-IR109

One-touch pop-up Ear Probe Covers



Accu-sensor
Technology



Ear Measure



Age Optional



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



One Second
Reading

Model No: FC-IR109

Product dimension: 35*69.4*145.5mm

Packing size: 45*80*165mm

Full set weight: 129.4g

Weight of thermometer: 82.4g(Without battery)/85.6g

Quantity per carton: 100 Pieces

N.W/Carton: 14.5kg

G.W/Carton: 16kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box



Infrared Thermometer

CE FDA



Package includes Ear
Thermometer Probe Covers

Ear Thermometer Probe Covers Test



FC-IR109

FC-IR109LED



Accu-sensor
Technology



Ear Temperature
Measure



Baby Mode



One Second
Reading



Three Color
Indicator Light



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



LED display

One-touch pop-up Ear Probe Covers



Model No: FC-IR109LED

Product dimension: 35*69.4*145.5mm

Packing size: 45*80*165mm

Full set weight: 129.4g

Weight of thermometer: 82.4g(Without battery)/85.6g

Quantity per carton: 100 Pieces

N.W/Carton: 14.5kg

G.W/Carton: 16kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box



Infrared Thermometer

CE FDA



FC-IR202

FC-IR202



Accu-sensor
Technology



Forehead
Measure



Light Spot Positioning



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



Model No: FC-IR202

Product dimension: 36*42*153.5mm

Packing size: 46*46*168mm

Full set weight: 115g

Weight: 66.8g(Without battery)/81.4g

Quantity per carton: 100 Pieces

N.W/Carton: 12.5kg

G.W/Carton: 14kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch,
2*Batteries, 1*Color box



Infrared Thermometer

CE FDA



FC-IR202LED

FC-IR202LED



Accu-sensor
Technology



Forehead
Measure



Baby Mode



One Second
Reading



Three Color
Indicator Light



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



LED display



Model No: FC-IR202LED

Product dimension: 36*42*153.5mm

Packing size: 46*46*168mm

Full set weight: 115g

Weight: 66.8g(Without battery)/81.4g

Quantity per carton: 100 Pieces

N.W/Carton: 12.5kg

G.W/Carton: 14kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch,
2*Batteries, 1*Color box



Infrared Thermometer

CE FDA



FC-IR1010LED

FC-IR1010LED



Accu-sensor
Technology



Ear/Forehead
Measure



Baby Mode



One Second
Reading



Three Color
Indicator Light



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



LED display

Model No: FC-IR1010LED

Product dimension: 36*42*153.5mm

Packing size: 46*46*168mm

Full set weight: 115g

Weight: 66.8g(Without battery)/81.4g

Quantity per carton: 100 Pieces

N.W/Carton: 12.5kg

G.W/Carton: 14kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch,
2*Batteries, 1*Color box

Infrared Thermometer


CE FDA



FC-IR206

FC-IR206

Tips for temperature measurement



Milk



Water



Room



Food



Accu-sensor
Technology



Forehead
Measure



Fever Alarm



Light Spot Positioning



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



One Second
Reading

Model No: FC-IR206

Product dimension: 150*38*50mm

Packing size: 46*46*168mm

Full set weight: 116g

Weight of thermometer: 70g(Without battery)/85g

Quantity per carton: 100 Pieces

N.W/Carton: 12kg

G.W/Carton: 14kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box

Infrared Thermometer

CE FDA



FC-IR207

FC-IR207



Accu-sensor
Technology



Forehead
Measure



Light Spot Positioning



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



Measure milk



Measure forehead



Measure room



Measure food

Model No: FC-IR207

Product dimension: 35*36*157mm

Packing size: 46*46*168mm

Full set weight: 136g

Weight: 69g(Without battery)/85g

Quantity per carton: 100 Pieces

N.W/Carton: 9kg

G.W/Carton: 10.5kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box

Infrared Thermometer

CE FDA



FC-IR207 V2

FC-IR207 V2



Accu-sensor
Technology



Forehead
Measure



Age Optional



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



Measure milk



Measure forehead



Measure room



Measure food

Model No: FC-IR207 V2

Product dimension: 35*36*157mm

Packing size: 46*46*168mm

Full set weight: 136g

Weight: 69g(Without battery)/85g

Quantity per carton: 100 Pieces

N.W/Carton: 9kg

G.W/Carton: 10.5kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box

Infrared Thermometer

CE FDA



FC-IR207LED

FC-IR207LED



Accu-sensor
Technology



Forehead
Measure



Baby Mode



One Second
Reading



Three Color
Indicator Light



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



LED display

Model No: FC-IR207LED

Product dimension: 35*36*157mm

Packing size: 46*46*168mm

Full set weight: 136g

Weight: 69g(Without battery)/85g

Quantity per carton: 100 Pieces

N.W/Carton: 9kg

G.W/Carton: 10.5kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Pouch
2*Batteries, 1*Color box



Infrared Thermometer

CE FDA



FC-IR2000

FC-IR2000



Accu-sensor
Technology



Forehead
Measure



One Second
Reading



°F/°C
Switchable



35 Groups Memory
Function



Object/Room
Measure



Fever Alarm



Tri-color Backlit

Model No: FC-IR2000

Product dimension: 105*35*137mm

Packing size: 116*41*150mm

Full set weight: 197g

Weight of thermometer: 79g(Without battery)/104g

Quantity per carton: 80 Pieces

N.W/Carton: 12kg

G.W/Carton: 13.5kg

Packaging: 1*Infrared Thermometer, 1*Manual, 1*Blister,
2*Batteries, 1*Color box



Finicare, a step towards health care.

Keep in Touch



+86 755 23281360



info@finicare.com



www.finicare.com



**201, No. 50, the 3rd Industrial Park Houting Community, Shajing Street
Bao'an District 518104 Shenzhen**